



# **Respiratory Monitoring**

**Ashraf Roshdy**

**Overview**

**Introduction**

**Clinical Monitoring**

**Lung Volumes and Pulmonary Function Tests**

**Monitoring Oxygenation / Ventilation**

**Imaging Monitoring**

**Monitoring during Mechanical Ventilation**

**Conclusion**

**References**

## **Introduction**

Respiratory failure, both acute (ARF) and chronic, is a life threatening global challenge. It is a dynamic complex syndrome originating from diverse pathologies; acute respiratory distress syndrome (ARDS), affecting a quarter of ventilated patients, being the most challenging. As unguided interventions were shown to add insult to the original injury, diagnosis and close monitoring are warranted. While remote community monitoring is growing, this chapter will focus on bedside respiratory monitoring in the acute settings.

Clinicians aim to treat acute respiratory failure (ARF) with minimal therapy associated harm. While society guidelines and lung protective strategies provide a framework of parameters, targets, and limits, their application to individual patients can yield variable outcomes. The heterogeneity of patient specific factors including host response, underlying pathology, and disease severity renders a unified approach less effective. Additionally, respiration involves a dynamic interplay of multiple physiological mechanisms and organ systems. Clinicians must diagnose the underlying cause, monitor deterioration to enable timely intervention, provide safe respiratory support when necessary, assess for improvements and readiness to wean the support. Throughout this process, an individualized approach is increasingly advocated, necessitating an effective and comprehensive monitoring strategy.

### **Monitoring as integral to respiratory support**

Acute respiratory support (RS) is used as a bridge for recovery or long term support. It can take place anywhere from pre-hospital to the intensive care units (ICUs) and includes conventional oxygen, high flow nasal oxygen (HFNO), mechanical ventilation (MV) and extracorporeal support. The application is not uniform worldwide, including doctors, nurses or respiratory therapists with variable degrees of competencies, especially out of hours.

The aim of RS is to improve the gas exchange and work of breathing while keeping to a minimum the harms arising from the support itself (e.g., dyssynchronies, lung and diaphragmatic injury). A complex decision making is needed in the more complex cases and the role of monitoring here is to support decisions to be precise, accurate, safe and timely.

In the ICU, precision is mostly physiologically based. Respiration involves the constant interplay of many organs (e.g., brain, heart, lungs, muscles, etc.) (Fig 1). Effective monitoring will capture, record, integrate, simplify and display such complexity and dynamism for every single patient. One step further is the move by the industry to integrate all monitored data into a single platform and furthermore use machine learning

and artificial intelligence to create autonomous ventilator adjustments in a closed loop fashion. Nevertheless, this remains currently premature, and its blind application is not without risks.

### Principles of Monitoring

Monitoring can be basic or advanced, intermittent or continuous. Targeting the gas exchange at its end point (e.g., blood or tissue oxygenation) is common but overlooks the underlying pathophysiology. In many cases, the therapy itself can influence many elements simultaneously, whether in a positive or a negative way. Identifying the underlying derangement(s), whether the respiratory drive, neuromuscular system, airways, lungs and chest wall mechanics, is the goal of advanced monitoring.

Monitoring can be clinical or requires sophisticated equipment. Mechanical ventilators themselves are increasingly serving as monitors but with some limits. An ideal monitor will provide a bedside, real time, continuous, accurate, inexpensive, simple to interpret, free from noise signals. Monitoring per se is not expected to improve the outcome, while the correct timely action by the operator can. As such, without a proper understanding and application, an expensive equipment can become useless or even harmful.




Figure 1: The different targets of respiratory monitoring

## Respiratory monitored parameters

### Clinical respiratory monitoring

#### 1. Respiratory drive

The respiratory drive refers to the intensity of the respiratory centers output (RCO) which mainly controls the ventilation (CO<sub>2</sub> elimination). RCO is then translated into an inspiratory flow-generation pathway through a neural pathway and the respiratory muscles. A depressed drive will lead to respiratory failure, while an increased one can lead to ventilator induced lung injury (VILI), self-inflicted lung injury (SILI) and ventilator induced diaphragmatic dysfunction (VIDD). Nevertheless, it is impossible to directly measure the respiratory drive, and it is usually assessed through its mechanical output (breathing effort).

The central control of breathing resides in the brain stem and cerebral cortex. The central control is connected to a control arm (flow-generation neuromuscular pathway) and feedback sensory mechanisms. Carbon dioxide is the strongest stimulus of the respiratory centers, however, the slope of the response can be affected by hypoxia, the acid base status, sleep and sedation.

During critical illness, it is possible that the RCO to be discrepant from the actual Minute Ventilation ( $V_{\min}$ ) if the control arm is dysfunctional or the lung mechanics are deranged. Vaporidi et al suggested a brain curve and ventilation curve to illustrate such discrepancy. The 2 curves are identical in physiological states, but when the ventilation curve shifts downward and to the right this leads to equilibrium at higher PaCO<sub>2</sub> levels. This discrepancy can be due to the respiratory center/brain curve (e.g., hyperoxemia, metabolic alkalosis, sleep and sedation) or inspiratory flow-generation pathway dysfunction (i.e., ventilation curve). To note that the respiratory drive can also recruit the expiratory muscles during respiratory distress. Apart from respiratory failure, assessing the respiratory drive is essential also during weaning from mechanical ventilation.

As mentioned above, the respiratory drive cannot be estimated directly, either the P<sub>0.1</sub> or the breathing effort (neuromuscular function) are used as a surrogate (see later).

## **2. Respiratory rate (RR)**

Patients recruit their respiratory reserve either by increasing the inspired volume (tidal volume) or by increasing the respiratory rate (RR).

RR is one of the four vital signs and a cornerstone of the early warning scores. Concomitant pathology (e.g., sepsis, acute neurological injury) and drugs (e.g., opiates) can interfere with the respiratory centers and lead to inappropriately increased or reduced RR. Despite such caveats, RR can precede changes in other early warning components and its extreme values are associated with mortality. In critically ill patients, the unstressed RR was found to be higher than normal adults by about 10 cycles/min. This is true also in patients with brain injury when they lose the cortical inhibitory effect exerted on the respiratory centers.

RR can be simply counted or watched on bedside monitors. More importantly, it is to recognize the underlying cause (whether respiratory decompensation or other causes), but unfortunately it was shown that clinicians are less mindful of the RR significance. The increase in spontaneous RR in respiratory failure most often follows the increase in muscular effort, so considered a sign of fatigue and decompensation which deserves prompt action.

In mechanically ventilated patients, the RR is set by the operator in controlled modes or triggered by the patient in assisted and spontaneous modes. Displayed RR might be different than the intrinsic patient's trigger. Ineffective efforts and false triggering are examples. Tachypnea can be a sign of under assistance. While the harm of pressure and volume are much appreciated by physicians, the contribution of a high RR (and flow) are generally less appreciated (see mechanical power). This includes hyperinflation (intrinsic PEEP) due to inadequate expiratory time, increase in the right ventricular afterload, and multiplying the times an injurious breath and energy are applied to the lungs (mechanical power).

## **3. Respiratory pattern and accessory muscle use**

In response to an increased respiratory drive, more forceful contraction of the muscles is needed to increase the inspired volume. The main muscle is the diaphragm, but there might be a need to recruit more muscles (accessory muscles). The excessive use of accessory muscles and expiratory muscles is a sign of deterioration. By increasing lung strain and stress, such excessive effort can lead to patient SILI and myotrauma.

## **Lung Volumes Monitoring**

### **Pulmonary Function Test (PFT)**

(Will be covered in detail in separate chapter)

In spontaneous breathing patients, Pulmonary function tests (PFTs) are crucial diagnostic tools used to assess the functionality of the respiratory system. These tests encompass a variety of measurements that evaluate how well the lungs take in and release air, and how efficiently they exchange gases and pulmonary impairment. PFT has diagnostic and therapeutic roles. Key parameters such as lung volume, capacity, airflow, and gas exchange efficiency are assessed through PFTs. Analyzing these metrics, aid in diagnose and monitor conditions like obstructive and restrictive lung diseases. PFTs play a pivotal role in guiding treatment decisions and tracking the progression of respiratory diseases.

### **Tidal volume (VT)**

During mechanical ventilation, Tidal volume (VT) is readily displayed on the invasive and non-invasive ventilators' screens. However, tidal volumes are not directly measured by the ventilator but integrated from the flow which is directly measured ( $VT = \text{Flow} \times \text{Time}$ ) usually the expiratory one (VTe). The difference between the inspiratory and expiratory tidal volume reflects the leak and sometimes the ventilator tubing compliance.

Since the ARMA trial, a VT of 4-8 ml/kg PBW (predicted body weight) was recommended in invasive MV in ARDS. VT causes lung expansion (deformation) which represent the lung strain when correlated to the baseline size (EELV or functional residual capacity, FRC). Even a slightly elevated VT can lead to volutrauma in the less aerated derecruited lungs (baby lungs). VT is difficult to control during non-invasive ventilation (NIV), and a value  $> 9.5$  ml/Kg PBW can predict NIV failure in hypoxemic respiratory failure.

### **Gas Exchange Monitoring**

The main function of the respiratory system is to secure an adequate gas exchange (i.e., provide adequate oxygen and eliminate carbon dioxide) to fulfil the metabolic needs. While lungs oxygenate blood passing through, carrying and extracting oxygen at tissue level requires adequate circulatory and cellular functions. As such, respiratory monitoring should not be a standalone approach when facing multi-organ dysfunction.

## 1. Pulse oximetry

Different hemoglobin (Hb) types absorb light at different wavelength. Pulse oximeters translate such phenomenon into a percentage of oxygen saturation, SpO<sub>2</sub> (OxyHb/total Hb). It captures also the pulse rate based on plethysmography, which when visually displayed reflects the quality of the signal and hence accuracy of the reading. Oximeters are widely used due to their simplicity and low cost and can be placed on the fingers, toes, forehead, nose or ears. Home oximeters are available, but with variable accuracies.

Skin color, perfusion, temperature and motion can impede the signal and lead to false results. There is a weak correlation with partial arterial oxygen tension (PaO<sub>2</sub>) at extreme values (e.g., SpO<sub>2</sub> <80% or 100%). The presence of abnormal Hb (e.g., Carboxyhemoglobin and Methemoglobin) can lead to false readings. To note that the numerical target of SpO<sub>2</sub> was recently revisited in many trials and best to be personalized with growing emphasis on the harm of hyperoxemia.

The SpO<sub>2</sub>/FiO<sub>2</sub> ratio was suggested as a surrogate of the classic PaO<sub>2</sub>/FiO<sub>2</sub>. This is of particular interest in resource limited settings or in out of ICU settings where less staffing is available and arterial lines not in use. The new global definition of ARDS supported using SpO<sub>2</sub>/FiO<sub>2</sub> ratio ≤ 315 for diagnosis (equivalent to PaO<sub>2</sub>/FiO<sub>2</sub> <300). In such way, SpO<sub>2</sub>/FiO<sub>2</sub> can more rapidly and accurately detect and classify ARDS (SpO<sub>2</sub>/FiO<sub>2</sub> thresholds: mild 234-314; moderate 148-235, or severe ARDS ≤ 148). Nevertheless, pulse oximeters are late detectors of hypoventilation (type 2 respiratory failure) and cannot stand alone as a respiratory monitor.

Lastly, monitoring the combination of RR and SpO<sub>2</sub> can predict the failure of HFNO. This can be achieved by using the ROX index:

$$\text{ROX Index} = (\text{SpO}_2/\text{FiO}_2) / \text{RR}$$

A lower value denotes more probability of HFNO failure, and a value greater than 4.88 was associated with less risk of progression to intubation. Different time points (ranging 2 and 12 hours after HFNO initiation) were used in different studies, with 6 hours point being the most common, showing low sensitivity but acceptable specificity (0.67 and 0.72 respectively).

## 2. Tissue oxygenation

While oximeters provide data about the peripheral arterial oxygen supply, discrepancies at tissue, cellular and mitochondrial levels do occur. Tissue oxygen tension (tPO<sub>2</sub>) refers to the oxygen partial pressure within the interstitial space of a particular tissue reflecting the balance between supply and utilization. Organ specific tools are available including the brain (PBtO<sub>2</sub>). Intuitively, the management should target perfusion, Hb concentration as well as the respiration. However, in many instances, the first, fastest and easiest response is to increase the fraction of inspired O<sub>2</sub> (FiO<sub>2</sub>) and hence pragmatically integrated as a target for any respiratory support.

Monitoring venous oxygen saturation (PVO<sub>2</sub>) is another way to assess the adequacy of oxygen supply and extraction. A lower concentration signifies decreased oxygen delivery or higher extraction which can be organ based (e.g., jugular venous oxygen saturation) or systematic (e.g., mixed venous oxygen saturation). Early during shock states, a mixed venous saturation of 70-80% is usually targeted despite downplayed by recent guidelines. Continuous monitoring is possible through fiberoptic catheter that uses reflection spectrophotometry.

## 3. Capnography

The second function of the respiratory system is to eliminate carbon dioxide (CO<sub>2</sub>). This can be measured in the peripheral blood (arterial blood gas) or in the exhaled CO<sub>2</sub> (PECO<sub>2</sub>) either in intubated or non-intubated patients. Integration into non-invasive devices also does exist.

Expired CO<sub>2</sub> is the endpoint of its production (metabolism), transport (circulation) and elimination (respiration). The arterial partial pressure of CO<sub>2</sub> (PaCO<sub>2</sub>) in the arterial blood gases (ABG) reflects the relation between CO<sub>2</sub> production and alveolar ventilation. Alveolar ventilation equals minute ventilation (V<sub>min</sub>) minus the dead space and represents the volume contributing to gas exchange. The dead space is the fraction of the tidal volume (V<sub>T</sub>) not participating in the gas exchange whether in the airways (anatomical), alveoli (alveolar) both representing the physiologic dead space in addition to the ventilator circuit (mechanical). The normal anatomic dead space equals approximately 2.2 ml/Kg.

$$\text{PaCO}_2 = 0.863 \times [ V_{\text{CO}_2} / V_A ]$$

Where V<sub>A</sub> is alveolar ventilation, V<sub>CO<sub>2</sub></sub> is CO<sub>2</sub> production, PaCO<sub>2</sub> is arterial partial pressure of carbon dioxide



$$V_A = V_{\min} \times (1 - V_d/V_T)$$

Where  $V_{\min}$  is minute ventilation,  $V_d$  is the physiologic dead space and  $V_T$  is the tidal volume

In a resting healthy adult,  $V_{CO_2}$  is approximately 200 ml/min and  $V_D/V_t$  is 0.3, which can be altered in the disease state.

The expired end tidal value ( $ET_{CO_2}$ ) is usually used to monitor the trend of the  $Pa_{CO_2}$ , while the difference between the 2 values can indicate the (physiologic) dead space. This is represented by the Enghoff's Modification of the Bohr Equation

$$\text{Alveolar Dead space}/V_T = (Pa_{CO_2} - P_{ET_{CO_2}}) / Pa_{CO_2}$$

Where  $P_{ET_{CO_2}}$  is the end tidal  $CO_2$ ,  $Pa_{CO_2}$  is the arterial partial pressure of  $CO_2$ ,  $V_T$  is the tidal volume

The dead space is an independent predictor of mortality in ARDS, and in order to avoid alveolar overdistention became a target for personalized mechanical ventilation. Otherwise, a trace (capnogram) can be displayed which helps to confirm the placement of endotracheal tube (ETT), diagnose bronchospasm, rebreathing and asynchrony, to mention few.




Figure 2: Schematic drawing shows 3 different alveoli with their corresponding capillaries. On the left: well perfused and ventilated alveolus with effective alveolar volume ( $VE_{\text{alv}}$ ), the middle is a ventilated but not perfused alveolus and is considered as alveolar dead space ( $VD_{\text{alv}}$ ), and the right one is an alveolus that is perfused but filled with fluid or collapsed (shunt). From reference 29.




Figure 3: Estimation of dead space and alveolar tidal volume from the volumetric capnometry waveform. The black dashed line a–b defines equal area q and p in phase II. The volume to the left of this line represents the anatomical dead space, while the right of the line is the alveolar tidal volume. The blue dashed line c–d is created so that area A equals area B. The distance from b to d defines alveolar dead space.  $V_{Daw}$ : anatomical dead space,  $V_{T_{alv}}$ : alveolar tidal volume,  $V_{D_{alv}}$ : alveolar dead space,  $V_{D_{phy}}$ : physiological dead space,  $V_{T_{alv-eff}}$  ( $V_{E_{alv}}$ ): effective alveolar volume,  $V_T$ : total tidal volume. Reference: Daoud EG, Franck CL. Alveolar mechanics: A New Concept in Respiratory Monitoring. From reference 29.

## **Role of Imaging in respiratory monitoring**

(Will be covered in detail in separate chapter)

With the exception of Electrical impedance tomography (EIT), imaging lacks the real time continuity of other monitoring devices. Moreover, computed tomography (CT) cannot be applied bedside and needs to carry the risk and workload of transfer. However, CT played a major role to delineate the pathophysiology and phenotypes of ARDS (e.g., baby lung) and still serve as a reference method to validate bedside recruitment tests. Moreover, lung and diaphragmatic ultrasound use is booming as a friendly bedside tool. Respiratory imaging will be revisited in more details in a dedicated chapter of this textbook.

### **1. Chest X-ray**

Chest X-ray (CXR) imaging is a widely used simple diagnostic tool for monitoring lung health and detecting various pulmonary conditions. It provides a two-dimensional image of the chest, allowing visualization of the lungs, heart, mediastinum, airways, blood vessels, and bones. This non-invasive technique is particularly valuable for identifying abnormalities such as pneumonia, lung tumors, tuberculosis, pleural effusion, and chronic obstructive pulmonary disease. While CXR are quick, cost-effective, and involve relatively low radiation exposure, they may not always provide detailed information about subtle or early-stage lung changes. In many instances also, the image quality can be an issue. In such cases, advanced imaging modalities like CT scans may be recommended. Despite its limitations, CXR imaging remains a cornerstone first line tool in lung monitoring, offering critical insights for diagnosis, treatment planning, and follow-up care.

### **2. Computed tomography (CT)**

CT of the chest is a highly advanced imaging modality that provides detailed cross-sectional images of the respiratory system, including the lungs, trachea, bronchi, and pleura. Utilizing X-rays and computer processing, CT chest imaging offers superior spatial resolution and contrast compared to conventional radiography, enabling the detection of subtle abnormalities such as small nodules, early interstitial lung disease, or vascular anomalies (e.g., thromboembolism). High-resolution CT (HRCT) is particularly valuable for evaluating diffuse lung diseases, as it can delineate fine parenchymal structures and patterns of pathology, such as ground-glass opacities, fibrosis, or emphysema. Contrast-enhanced CT is often employed to assess vascular structures, mediastinal masses, or infectious processes, while low-dose CT has become a cornerstone in lung cancer screening for high-risk individuals. The versatility and precision of CT chest imaging make it an indispensable tool in the

diagnosis, staging, and management of a wide range of respiratory conditions, from infections and malignancies to chronic obstructive pulmonary disease and pulmonary embolism. Moreover, the nature and distribution of the structural injury can guide interventions (e.g., proning and recruitment in ARDS, pleural drainage, bronchoscopy in lung collapse, etc.).

### **3. Magnetic resonance imaging (MRI)**

MRI of the respiratory system is a non-invasive imaging technique that utilizes strong magnetic fields and radio waves to generate detailed images of thoracic structures, including the lungs, mediastinum, and chest wall. Unlike CT, MRI does not involve ionizing radiation, making it a safer option for certain patient populations, such as pregnant women or those requiring repeated imaging. MRI excels in soft tissue contrast, allowing for superior visualization of mediastinal and chest wall pathologies, including tumors, vascular anomalies, and inflammatory conditions. Techniques such as diffusion-weighted imaging (DWI) and dynamic contrast-enhanced MRI (DCE-MRI) provide functional insights, aiding in the characterization of lesions and assessment of tumor perfusion. However, MRI has limitations in evaluating the lung parenchyma due to low proton density and susceptibility artifacts from air-tissue interfaces, which can reduce image quality. Despite this, advances in ultrashort echo time (UTE) sequences and hyperpolarized gas MRI are expanding its utility in assessing lung structure and function, particularly in diseases like cystic fibrosis or chronic obstructive pulmonary disease (COPD). While not a first-line modality for routine respiratory imaging, MRI plays a complementary role in specific clinical scenarios, offering unique diagnostic information without radiation exposure.

### **4. Ultrasound**

Ultrasound imaging of the respiratory system is a versatile, non-invasive, and radiation-free modality that is particularly useful for evaluating superficial thoracic structures, pleural pathologies, and guiding interventional procedures. The Ultrasound waves (US) become absorbed as they pass through the air; a major barrier for its use. Nevertheless, in lung collapse/consolidation, US becomes an excellent visualizing tool of the de-aerated lungs.

Additionally, US had shown superiority in detecting pleural diseases like pleural effusions, and pneumothorax compared to conventional x-ray imaging. The growing use, skills and availability of US makes it an attractive tool. Serial use can substitute its non-continuity compared to other monitors.

To note that many protocols integrate cardiac and lung ultrasound. Obviously, screening for pulmonary hypertension and right heart strain can guide the respiratory management in many cases (e.g., PEEP setting). Also, the contribution of left heart pathologies in respiratory failure should not be overlooked. Last, heart-lung interaction is a cornerstone in the assessment for fluid responsiveness during shock state commonly associated with acute respiratory failure.

### **5. Ventilation / Perfusion (V/Q) Scan**

V/Q scan is a nuclear medicine imaging technique used to evaluate lung function by assessing airflow (ventilation) and blood flow (perfusion) within the lungs. This test is particularly valuable in diagnosing pulmonary embolism (PE), as it can identify mismatched defects where areas of the lung are ventilated but not perfused, a hallmark of PE. The use of V/Q scans have dropped significantly with the advancement of CT scan technology.

### **6. Electrical impedance tomography (EIT)**

EIT is a non-invasive, radiation-free imaging technique that provides real-time monitoring of lung function by measuring changes in electrical impedance within the thoracic cavity. Different tissues have different electric conductivity and resistance (impedance). Meanwhile, air is a poor conductor. A belt and 16-32 electrodes attached across the 5th-6th intercostal transverse plane measure the surface voltage following the application of an alternating electric current. The captured impedance of a lung slice is subsequently translated into an image with different color scales on a monitor. The display shows a real time image of the lungs but allows also a view of the different zones to detail regional ventilation, collapse and overdistention. EIT has the advantage of being bedside and allowing real time regional imaging offering clinicians valuable guidance.

## Advanced monitoring during mechanical ventilation

### 1. Lung mechanics

Understanding lung physiology and mechanics is fundamental to properly set any respiratory support. The lungs are cyclically exposed to a change in volume. The volume change compared to baseline is called strain. The driving force behind the flow of air from the atmosphere into the lungs (ventilation) is a pressure gradient (driving pressure,  $\Delta P$ ). The lungs behave as an elastic organ which tend to collapse when expanded due to its elastance. Elastance is the change of pressure divided by the change of volume. Compliance (C) is the reciprocal of elastance ( $E=1/C$ ). The pressure represents lung stress which correlates to the strain, and both are implicated in ventilator induced lung injury (VILI) if certain limits are exceeded

$$\text{Stress} = K \times \text{Strain}$$

Where K is a constant (specific elastance)

The pressure drop across the respiratory system sums the drop during the air passage in the airways (between the airway opening and the alveoli) and the pressure needed to expand the lungs and the chest wall. An additional minimum pressure is needed to overcome the gas and tissue inertia (temporal acceleration of gas).

$$\Delta P = P_{\text{resistance}} + P_{\text{elastance}} + P_{\text{inertance}}$$

Inertia is usually negligible and omitted for simplicity.  $\Delta P$  represents the gradient between the pressure at the airway opening ( $P_{aw}$ ) and any end expiratory pressure within the alveoli (i.e., total PEEP).  $P_{\text{resistance}}$  can be calculated as the product of the gas flow multiplied by the airway resistance, while the elastance pressure is the product of the volume multiplied by the respiratory system elastance ( $E_{rs}$ ). All this is represented by the equation of motion:

$$P_{aw} = (R \times F) + (V \times E_{rs}) + PEEP$$

Or

$$P_{aw} = (R \times F) + (V/C_{rs}) + PEEP$$

Where R is resistance, F is flow,  $E_{rs}$  is respiratory system elastance,  $C_{rs}$  is respiratory system compliance and V is the volume.

The pressure can be generated by the patient's muscular effort, MV or a combination of both. As they are arranged in series, the generated pressures equal their sum. In the spontaneously breathing patient, the muscular contraction generates a negative intrapleural pressure sucking the air from the atmosphere (negative pressure ventilation). In controlled positive pressure ventilation, the ventilator creates a positive pressure at the airway opening (external port of the endotracheal tube) and the gas is pushed to the alveoli. In assisted mechanical ventilation the patient effort adds to the pressure generated by the ventilator.

$$P_{aw} = P_{mus} + P_{vent} = (R \times F) + (V/C_{rs}) + PEEP$$

Where R is resistance, F is flow, C is respiratory system compliance, and V is volume.

In static conditions of no flow, the contribution of resistance is abolished from the equation (as flow equals zero). This is the reason many indices are measured after applying an inspiratory and expiratory hold (static measurements). When performing an inspiratory hold, the pressure drop from a peak to a plateau pressure ( $P_{plateau}$ ) reflects the  $P_{resistance}$ . It becomes then possible to calculate the  $P_{elastance}$  by subtracting the PEEP from the  $P_{plateau}$  (yielding the driving pressure;  $\Delta P$ ). Compliance (static) can then be calculated by dividing VT by  $\Delta P$ .

$$C_{rs} = VT / \Delta P$$

$\Delta P$  reflects the normalized tidal volume to the respiratory system compliance ( $C_{rs}$ ); the later in turn correlates with the size of aerated lungs and expectedly associated with ARDS mortality even when delivering  $V_T$  or  $P_{plateau}$  within the protective range.  $\Delta P$  is generally recommended to be kept  $< 15$  cmH<sub>2</sub>O, and one work showed a 5% increase in mortality with every subsequent increment of 1 cmH<sub>2</sub>O.

It was recognized for long time that a high pressure applied to the alveoli can cause injury and air leak (barotrauma), and hence the recommendation to limit  $P_{plateau}$ . Dreyfus et al subsequently showed that a high  $P_{aw}$  is less injurious if achieved by strapping the chest. A significant proportion of  $P_{aw}$  in this case was from the external pressure exerted by strapping leading to high pleural pressure ( $P_{pl}$ ). This led to the better understanding that only the pressure needed to distend the lung tissue (i.e., the transpulmonary pressure (TPP) is responsible for the lung stress and VILI rather than the total pressure applied to the whole respiratory system: lungs and chest wall (i.e.,  $P_{aw}$ ).

$$\mathbf{TPP = P_{alv} - P_{pl}}$$

$P_{pl}$  is dependent on the chest wall elastance ( $E_{cw}$ ). The total respiratory system elastance sums the lung elastance ( $E_L$ ) and the chest wall elastance ( $E_{cw}$ ).

$$\mathbf{E_{rs} = E_L + E_{cw}}$$

TPP cannot be simply estimated from  $P_{aw}$  as the lung elastance varies and can range anywhere between 50 to 90% of the total  $E_{rs}$ . For example, in ARDS, the pressure needed to oppose the chest wall elastance is in the range of 5–10 cmH<sub>2</sub>O but can vary significantly in obese patients or in case of high intra-abdominal pressure. To this end, it is clear that partitioning the distending pressure of the lungs from the chest wall is of major clinical interest. To achieve that, the simplest way is to directly measure the pleural pressure:

$$\mathbf{P_{cw} = P_{pl} - P_{bs}}$$

Where  $P_{cw}$  is the chest wall distending pressure,  $P_{pl}$  is pleural pressure and  $P_{bs}$  is the pressure at the body surface.

While this is not feasible in clinical practice, the pressure in the adjacent esophagus ( $P_{es}$ ) can be measured by the mean of an esophageal balloon inserted in its lower third (around 35-45 cm from the nostrils) and was widely validated as a surrogate of the  $P_{pl}$ . TPP can then be calculated as the difference between the airway pressure and  $P_{pl}$  indicating the lung parenchymal stress responsible for lung injury induced either by the ventilator (VILI) or the patient himself (SILI).

$$\mathbf{TPP = P_{alv} - P_{pl}}$$

Or

$$\mathbf{TPP = P_{alv} - P_{es}}$$






Figure 4: Left graph from a volume-controlled mode: Airway pressure (top in yellow) end expiratory airway pressure 15, peak inspiratory pressure 28, plateau inspiratory pressure 25. Flow (2nd row in pink): continuous flow of 20, end inspiratory and end expiratory flow of zero. Esophageal pressure (3rd row in orange): end expiratory pressure 14, peak inspiratory pressure 16, end inspiratory pressure 14. Airway pressure – Esophageal pressure (bottom row in orange) trans-pulmonary pressure of 12, end inspiratory trans-alveolar pressure 11, end expiratory trans-alveolar pressure 1. Right chart: V<sub>Daw</sub>: anatomical dead space, V<sub>Talv</sub>: alveolar dead space, Dead space/Tidal volume: V<sub>Daw</sub>/VTE Pressures in cmH<sub>2</sub>O, Flow in L/min, Volumes in ml. From reference 58.

## 2. Esophageal manometry

The principle for Esophageal manometry is described above. Esophageal manometry can serve to:

- i. Set the ventilator settings
- ii. Monitor the work of breathing (WOB)
- iii. Detect patient-ventilator interactions and asynchrony
- iv. Research tool to validate less invasive bedside tests

The partition of lung mechanics (pressure and compliance) and measuring TPP are recommended to personalize MV settings and avoid VILI especially in severe ARDS. In such cases, P<sub>es</sub> is used to personalize P<sub>aw</sub> and PEEP.

The pleural pressure is not uniform and can vary in the antero-posterior axis, and the directly measured  $P_{es}$  was found to correlate more with the pressure in the dependent lung zones. The elastance method to calculate  $P_{es}$  can correlate better with the non-dependent areas of the lung.

Inspiratory TPP remains a positive value whether during MV or spontaneous breathing. In the former MV generates a positive  $P_{aw}$  while in the later the patient's muscular effort generates a negative  $P_{pl}$ . In assisted MV modes, TPP is the sum of both. Vigorous muscular contraction increases the swing the  $P_{es}$  (reflecting the muscular pressure  $P_{mus}$ ) and a value  $\geq 15$  cmH<sub>2</sub>O can be used as a marker of increased work of breathing (WoB) and risks SILI. This can guide the escalation from non-invasive respiratory support including intubation.

The goal during inspiration is to avoid excessive stress (TPP), while during expiration is to avoid atelectasis and atelectotrauma by avoiding negative values.

$$\mathbf{TPP_{end - inspiratory} = P_{plat} - P_{es}}$$

$$\mathbf{TPP_{end - expiratory} = PEEP - P_{es}}$$

(measurements during end-inspiratory and end-expiratory pause respectively)

There is no definite consensus about  $TPP_{end-inspiratory}$ , but values  $\leq 15-20$  cmH<sub>2</sub>O can be reasonable. As mentioned above, the difference between the end inspiratory and end expiratory pressure ( $P_{plateau} - PEEP$ ) is called the driving pressure ( $\Delta P_{rs}$ ) and was correlated with mortality in ARDS. It combines the effect of the tidal volume and compliance/elastance.

$$\mathbf{\Delta P = VT \times E = VT/C_{stat}}$$

Measuring the driving TPP ( $\Delta P_L$ ) can more accurately reflect the tidal lung stress and guide MV. A safe limit in the range 10-12 cmH<sub>2</sub>O was suggested by experts but with little consensus.

Unfortunately, and despite its fundamental role, esophageal pressure monitoring is applied in only about 1% of ARDS patients.

### **3. Ventilator parameters**

Modern ventilators may be the most widely used respiratory monitors. In fact, manufacturers compete to offer advanced monitoring options and integrate them into algorithms for ventilator self-actions.

Essentially, the MV operator sets input variables which are translated (through lung mechanics) into output variables (uncontrolled by the operator). Modes are generally divided into volume and pressure modes considered as 2 faces of the same coin, where the relation between the volume and the pressure is governed by lung mechanics (equation of motion).

The displayed data can be set or measured parameters (e.g., RR, VT and pressures) with further calculation of the compliance, resistance and sometimes more advanced parameters (e.g.,  $P_{0.1}$ , leak and intrinsic PEEP). The continuous visual display of the volume, pressure and flow against time as curves (scalars) and against each other: the pressure-volume and flow-pressure loops, are becoming a standard in many ventilators.

The overall aim is to guide the operator to adjust the settings to match the patient's needs, avoiding harm and patient-ventilator asynchrony (PVA). Intelligent closed loop modes are emerging where the ventilators adjust the settings automatically based on the monitored data through pre-programmed algorithms, and the role of clinicians is to set the overall strategy then supervise its application.

Flow curves and loops are particularly useful to detect any mismatch or asynchrony between the patient and ventilator.

#### **4. Alveolar opening and closure pressures**

At static states (no flow), the  $P_{aw}$  and PEEP reflect the alveolar pressure only if there is a continuity between the proximal airway and the alveoli. In fact, it was shown that around a fifth of ARDS and obese patients exhibit an airway opening pressure (AOP) delineating small airway collapse during expiration. This highlights that collapse can happen not only at the alveolar but also the bronchiolar level leading to a discrepancy between the lungs macro-parameters, as targeted by MV settings, and the micro-environment at the alveolar level. This can lead to misleading measurements even with an expiratory occlusion pressure.

Different alveoli have different opening and closing pressure. The critical closing pressure is the pressure at which alveoli will collapse, while the critical opening pressure (AOP) is the one at which it will (re)open; the last being higher than the former (AOP > closing pressure). In order to avoid atelectotrauma secondary to alveolar cyclic opening and closing, PEEP should be set to keep alveoli open. Ideally, this can be achieved by setting PEEP higher than the closing pressure, despite in clinical practice the AOP is usually the used target, despite some risks of overdistention.

AOP can be detected on the pressure time curves or pressure-volume loops. Using volume mode and slow respiratory rate, the sudden change in the slope of the inspiratory ascending limb (pressure-time curve) signifies the sudden change of compliance from ventilator circuit compliance to the Crs as the alveoli open when reaching the AOP. On a pressure-volume loop, AOP corresponds to the sharp bend of the inspiratory limb (Figure 5). To note, that as recruitment is time as well as pressure dependent, the opening pressure can be better represented by a range of pressure rather than a single value.




Figure 5: Quasi-static Pressure-Volume curve (top green). Hysteresis volume (bottom yellow). Red arrow points to airway opening pressure, red circle point to point of maximum curvature.

## 5. Recruitability

‘Baby lung’ is a term used to describe the lung heterogeneity in ARDS where some pulmonary regions maintain normal inflation while the majority are collapsed. This led to an increased interest in recruiting the collapsed lung units to serve 2 goals:

- (1) Improving the gas exchange: more alveoli participating in gas exchange will reduce shunting through the lung, a main mechanism of hypoxemia.
- (2) Reducing VILI: When the VT distributes across more units (alveoli), it will produce less strain of each of them. Obviously, keeping more alveoli aerated at the end of expiration means an increase in the end expiratory lung volume (EELV), the denominator of the strain equation (see above).

To achieve that, alveoli need to be opened (inspiratory pressure exceeding the AOP), then prevented from collapsing during expiration (PEEP above the closing pressure). Some alveoli are slower in opening and closing than others (*see time constant*), and as some alveoli are already open and the intervention is applied to the lung as a whole, a recruitment maneuver (RM) can end with either recruitment or overinflation. A third scenario will take place in the consolidated regions (alveoli flooded with exudate or fluid) which will remain closed (resistant to recruitment). If the alveoli then re-close during expiration (PEEP below closing pressure), this will lead to injury by cyclic opening and closing (atelectotrauma).

RM responsiveness is highly variable and difficult to predict by static parameters. One work showed the range of recruitment to be 0 to 40% when increasing the PEEP from 5 to 10 cmH<sub>2</sub>O. More importantly, recruitment trials were negative and showed even an increased mortality in compliance based RM combined with high PEEP. Due to those reasons, RM universal application (without pretesting for recruitability) cannot be advised.

CT scan played a major role to describe the ‘baby lung’ and quantifying the pulmonary tissue and gas volume. The morphologic aeration was then correlated with compliance. While far from being a bedside test (apart from exceptional portable CT scanners), it led to a better understanding of ARDS pathophysiology and offered a tool to validate bedside tests.

### **6. Recruitment-to-Inflation ratio (R/I ratio)**

The concept is to calculate the (de)recruited volume ( $V_{rec}$ ) during a PEEP reduction (classically from 15 to 5 cmH<sub>2</sub>O). By using volume mode with slow respiratory rate (6-8 breaths/min), the  $VT_e$  following the PEEP drop will be larger than the previous  $VT_e$  as it sums the (de)recruited volume ( $V_{rec}$ ) (volume lost due additional alveolar collapse if derecruitment happens) plus the reduction in EELV at the lower PEEP (which can be expected by multiplying the  $\Delta P \times C_{rs}$  at the low PEEP).  $V_{rec}$  can then be calculated as the other variables are known ( $VT_e$ ,  $\Delta P$  and  $C_{rs}$ ) followed by calculating the compliance of the (de)recruited lung ( $V_{rec}/\Delta P_{PEEP}$  drop). R/I ratio is obtained by dividing the (de)recruited lung compliance by the  $C_{rs}$  at the lower PEEP. A value  $> 0.5$  was considered a marker of recruitability. While easy to perform bedside, the value can be over or underestimated using different ventilators, and correlate modestly with morphological recruitment on CT scans.

## 7. Hysteresis and Pressure volume (P/V) curve

Observing the pressure-volume loop, It is easy to notice that the inspiratory and expiratory limbs are not identical (Figure 6). The surface area between the two limbs is called hysteresis demonstrating that the pressure at a given volume is less during expiration than inspiration.

Hysteresis was correlated with compliance ( $C_{rs}$ ) and the dissipated energy within the lung tissue during respiration. It was postulated the larger the hysteresis the more the lungs' recruitability. This was supported by CT studies showing that hysteresis volume normalized to the maximum volume/pressure of the lung correlates with gas and tissue recruitment. It might be also that this dissipated energy contributes to molecular micro-injury (VILI), a point which needs further investigation.

Some ventilators introduced the P/V method (protective ventilation tool) (Figure 6). coupled by the option to perform a subsequent recruitment maneuver. Using low flow (constant pressure ramp 2-5  $\text{cmH}_2\text{O/s}$ ), the ventilator measures different parameters including the hysteresis volume difference at 20  $\text{cmH}_2\text{O}$  and the maximum one. The normalized maximal distance (NMD) is the maximum volume difference divided by the maximum volume ( $V_{\text{max}}$ ) reached at maximum pressure (usually 40  $\text{cmH}_2\text{O}$ ).  $\text{NMD} > 41\%$  indicates recruitability and the operator can proceed with RM through a pre-defined algorithm.

To note that both methods (RI and P/V) correlate well and are accepted bedside test for recruitability.




Figure 6: P/V loop showing Normalized Maximum distance and hysteresis. From Hamilton Medical with permission. In this case maximum  $dV=674$  ml, while maximum volume equals 1712 ml.  $\text{NMD} = 674/1712 = 39.4\%$ . From reference 80.

## **8. Positive End Expiratory Pressure (PEEP)**

The role of PEEP is to keep the alveoli open during expiration. Applying PEEP can be incremental or decremental (after RM) and needs to be set above the alveolar closing pressure, but the best way to identify the optimum PEEP remains under investigation.

Shortly after the description of ARDS, it was noticed that PEEP was associated with better oxygenation. This was attributed to alveolar recruitment, optimization of gas exchange, reduction in shunt fraction and lung strain (by increasing the end expiratory lung volume (EELV) the denominator of the strain equation =  $VT/EELV$ ), and the avoidance of tidal atelectotrauma (cyclical opening and closing of the alveoli). On the other hand, PEEP can lead to overdistention, increases the right heart afterload, and pushes the diaphragm down in a less optimum position for contraction. Such factors should be taken into account when setting the PEEP pointing to the complementary role of Critical Care Echocardiography in respiratory monitoring.

Optimum PEEP was a research topic for decades. One recent review identified 22 different investigated methods. The most common are the compliance-based, imaging-based (including EIT, CT and US) and PEEP-FiO<sub>2</sub> tables. Tables are the most simple, categorized mainly into high or low PEEP ones, but lack personalization. Another two common methods are a stepwise increment or decrement to reach the PEEP associated with best compliance or oxygenation (PaO<sub>2</sub> or SpO<sub>2</sub>).

In case of bronchospasm and/or short expiratory time, the alveoli do not empty completely leading to the presence of intra-alveolar air (and pressure) at the end of expiration. This is called intrinsic PEEP and is different than the PEEP set by the MV operator. It can be spotted if the expiratory flow does not reach the zero line at the end of expiration (flow-time curve) and can be quantified by performing an expiratory hold. Intrinsic PEEP imposes an additional WOB and can lead to asynchrony by making triggering more difficult (i.e., ineffective efforts). Bronchodilators, selecting the correct size of the endotracheal tube (ETT), and setting the correct expiratory time and inspiratory volume are needed to avoid intrinsic PEEP.

## **9. Time constant**

Time constant ( $\tau$ ) is the product of resistance and compliance. One time constant equals the time needed to inflate or deflate 63% of the total volume.  $3\tau$  and  $5\tau$  are needed to empty 95% and 99% of the volume and hence can be used to set the expiratory time.

$$\tau = C \times R$$

Where C is compliance (L/cmH<sub>2</sub>O) and R is resistance (cmH<sub>2</sub>O/L/s)

$\tau$  is about 0.1 seconds in a normal person. In pathological situations, an increase in compliance or resistance can increase  $\tau$ . PEEP also was shown to affect the  $\tau$ . As different alveoli can have different compliance or deflate against different airway resistance, they do have different time constant (i.e., they do not deflate at the same rate; high resistance for example will lead to prolonged time constant and slower filling and deflation).

### 10. Stress index (SI)

Observing the pressure-time curve during constant flow volume modes in a passive patient (no spontaneous effort) can help to assess lung recruitment or overdistention. An optimum inflation is represented by a linear slope of the inspiratory pressure-time curve. This can be assessed in a qualitative or quantitative way (stress index). An upward concavity is referred to as SI >1 and signifies overdistention. The opposite is a downward concavity (SI <1) indicates tidal recruitment. In such cases, PEEP may need to be adjusted. (Figure 7).




Figure 7: Airway Pressure-Time scalar using the volume controlled mode with constant flow in passive conditions showing different Stress Indexes. The red dashed line represent SI of 1.

### 11. Mechanical power

VILI involves the contribution of many variables. While the general recommendation is to limit the VT, P<sub>plateau</sub> and  $\Delta P$ , a global summative index including all those variables in one equation may be simply more informative. Energy is delivered to the lungs during breathing, but part of it dissipates within the lung tissue (see hysteresis) and is believed to be indulged in the molecular damage of the lung parenchyma.

Experts suggested the mechanical power (MP) as an inclusive physical measure to quantify the driving force for VILI. The equation was founded from the motion equation then multiplied by the RR to describe the energy delivered per minute (j/min). It integrates also the PEEP, but its contribution as



static pressure remains debated. The equation is not simple to calculate but provides a more global picture rather than the single variables. It can be calculated during invasive as well as NIV.

MP validity is based on many studies showing correlation with outcome in different ICU population including non-ARDS. The safe upper limit of MP was suggested to be 12-17 j/min at different time points in different studies. To note that adaptive ventilation modes were shown to be associated with less MP which points to a possible superiority in guarding against VILI. Some ventilators introduced a surrogate or simplified equations to calculate and display MP in volume and pressure modes.

## **12. Patient-ventilator asynchrony (PVA)**

Synchrony between the MV and the patient is essential not only to secure comfort but also for faster weaning and better survival. PVA arises from a mismatch between ventilator assistance and patient's needs either in the form of over or under assistance.

It needs expertise to detect, sometimes only possible with advanced monitoring (e.g., esophageal balloon). It can be simplified by considering the different phases of the respiratory cycle:

### **Trigger:**

- Delayed trigger
- Auto-triggering
- Ineffective effort

### **Flow:**

- Flow asynchrony (flow starvation)

### **Cycling (termination of breath):**

- Early cycling
- Delayed cycling
- Double triggering
- Reverse triggering

Ineffective effort is the most common followed by double triggering which can be more injurious due to the unintentional increase in VT (breath stacking). Identification is the first step to intervene, mostly by adjusting the ventilator settings. Deeper sedation can be sometimes the last resort but can lead or worsen reverse triggering.

# Respiratory Physiology – Respiratory Monitoring








Asynchrony	Description	On the waveform	Waveform example	Common possible causes
<b>Trigger asynchronies - during the beginning of inspiration</b>				
Delayed triggering	The time interval between the patient's inspiratory effort and the delivery of a mechanical breath is increased	Flow waveform: Look for a longer-than-normal time interval between the positive deflection in flow <b>1</b> and the delivery of ventilatory support <b>2</b>		<ul style="list-style-type: none"> <li>Trigger threshold set too high</li> <li>Ventilator pneumatics</li> <li>Presence of AutoPEEP</li> <li>Low respiratory drive</li> <li>Weak inspiratory effort</li> </ul>
Ineffective effort	The patient's inspiratory effort fails to trigger the delivery of a mechanical breath	Flow waveform: Look for an abrupt change in the steepness of the waveform <b>1</b> (decrease in expiratory flow or increase in inspiratory flow) that is not followed by ventilatory support <b>2</b>		<ul style="list-style-type: none"> <li>Trigger threshold set too high</li> <li>Pressure support too high</li> <li>Set frequency and/or inspiratory time too high (in controlled modes)</li> <li>Tidal volume set too high</li> <li>Presence of AutoPEEP</li> <li>Low respiratory drive</li> <li>Weak inspiratory effort</li> <li>Sedation</li> </ul>
Auto triggering	A mechanical breath delivered without an inspiratory effort	Pressure waveform: Look for a delivered mechanical breath showing no drop in airway pressure <b>1</b> at the beginning of the inspiratory phase		<ul style="list-style-type: none"> <li>Trigger threshold set too low</li> <li>Air leaks in the endotracheal tube cuff, ventilator circuit, or chest tube</li> <li>Flow oscillations (water or secretion in the circuit, cardiac oscillations)</li> </ul>
<b>Flow asynchronies - during the gas delivery</b>				
Flow asynchrony	The delivered flow does not meet the patient's inspiratory flow demands	Pressure waveform: Look for an upward concavity <b>1</b> preceding the end of the mechanical breath		<ul style="list-style-type: none"> <li>Inappropriate selection of ventilation mode (more frequent in volume-controlled modes)</li> <li>High inspiratory effort</li> <li>In volume-controlled modes: inappropriate flow settings</li> <li>In pressure-controlled modes: inappropriate P-ramp settings</li> </ul>
<b>Termination asynchronies - during the end of inspiration</b>				
Double triggering	Two (or more) mechanical breaths are delivered during one single inspiratory effort	Flow waveform: Look for two assisted breaths without expiration between them or with an expiration interval of less than half of the mean inspiratory time (often visually displayed as a waveform with two inspiratory peaks)		<ul style="list-style-type: none"> <li>Cycling criteria (ETS) set too high</li> <li>Pressure support too low</li> <li>P-ramp too short</li> <li>Flow starvation</li> <li>High respiratory drive</li> <li>Time constant too short</li> </ul> <p>Double triggering can be an effect of and/or promoted by reverse triggering or early cycling</p>
Early cycling	The duration of the mechanical breath is shorter than the duration of the patient's inspiratory effort	Flow waveform: Look for a small bump <b>1</b> at the beginning of expiration (after peak expiratory flow) followed by an abrupt initial reversal in the expiratory flow <b>2</b>		<ul style="list-style-type: none"> <li>In pressure support ventilation: <ul style="list-style-type: none"> <li>Cycling criteria (ETS) set too high</li> <li>Low levels of ventilator pressure support</li> </ul> </li> <li>Time constant too short</li> <li>In time-cycled ventilation: <ul style="list-style-type: none"> <li>Short inspiratory time</li> </ul> </li> </ul>
Delayed cycling	The duration of the mechanical breath is longer than the duration of the patient's inspiratory effort	Flow waveform: Look for a change in the slope of the inspiratory flow: a fast decrease <b>1</b> followed by an exponential (less steep) decline <b>2</b>		<ul style="list-style-type: none"> <li>In pressure support ventilation: <ul style="list-style-type: none"> <li>Cycling criteria (ETS) set too low</li> <li>Pressure support too high</li> <li>P-ramp too long</li> </ul> </li> <li>In pressure control ventilation: <ul style="list-style-type: none"> <li>Cycling criteria (ETS) set too low</li> <li>Inspiratory time too long</li> </ul> </li> <li>In volume control ventilation: <ul style="list-style-type: none"> <li>Low flow</li> <li>Long inspiratory time</li> <li>High tidal volume</li> </ul> </li> </ul>

Figure 8: Example of ventilator software recognition of dyssynchronies, their recognition, and causes. With permission and credit to Hamilton Medical, Bonaduz, Switzerland. From reference 64.

## 13. Leak

Leak is the difference between the delivered inspired VT ( $VT_i$ ) and expired VT ( $VT_e$ ). It is easily measured by the ventilators flow sensors and can happen at the ventilator circuit, endotracheal/tracheostomy tube or at the patient-ventilator interface in NIV. Sometimes, a false leak alarm is due to a dysfunctional sensor.

Leak is mostly inevitable during non-invasive ventilation (NIV) happening around the interface (or from the mouth in case of nasal mask). NIV can compensate for leaks, which is sometimes intentional to help clear CO<sub>2</sub>, but a large leak can end with less respiratory support (providing less pressure and VT), patient discomfort and interferes with triggering leading to PVA. This can end by failure of the NIV trial. While NIV has some capacity to compensate for leaks, this differs from device to device. The compensation is better in pressure targeted ventilators but can vary significantly among different devices. Eventually an intervention by the operator is required in the form of adjusting or replacing the masks or modifying the settings.

The endotracheal and tracheostomy tubes have a high-volume low-pressure cuff to guard against air leak and aspiration. During invasive MV, a leak should trigger first an assessment of the correct placement of the endotracheal or tracheostomy tube and the cuff pressure. Cuff pressure should be monitored (target: 20-30 cmH<sub>2</sub>O) every 8-12 hours to assure an effective seal but also to avoid pressure damage to the tracheal wall and interfering with its blood flow. New generation of cuff monitors reduce the nursing workload and continuously monitor and adjust the pressure (e.g., Intell cuff, cuff pressure manager). Those monitors can be standalone or integrated into the ventilator itself. A recent meta-analysis showed a benefit in terms of ventilator associated pneumonia (VAP), duration of MV and length of stay in the ICU.

Translaryngeal intubation can lead to laryngeal edema and subsequent failure of extubation in about 3.5% of patients. A cuff leak test is sometimes used to assess the risk of stridor post-extubation by demonstrating either audible leak or observing its volume.

### **14. Work of breathing**

The Work of breathing (WOB) was historically monitored as a sign of respiratory distress and the need to escalate support. WOB correlates with the energy expenditure of the respiratory muscles and is the product of the volume and pressure

$$\mathbf{WOB = \int \mathbf{Pressure} \times \mathbf{Volume}}$$

Paralleling the equation of motion, it sums 3 components: the work needed to overcome the (1) resistive force, (2) elastic force of the respiratory system and (3) intrinsic PEEP within the alveoli. It ranges between 2.4 to 7.5 J/min or 0.2 to 0.9 J/L . Some recent ventilators display such value based on the pressure time product (PTP) which reflects the pressure generated by the patient's inspiratory

muscles. A mismatch between the respiratory muscle load and capacity will end in respiratory failure unless assistance is applied.

A high central respiratory drive will manifest as an increased WOB. This can reflect the respiratory load (e.g., High  $\text{CO}_2$ , metabolic reasons) but sometimes it is unrelated (e.g., head injury, pain, anxiety, etc.). This will translate into an increase in the minute volume, through an increase in the muscular effort (and hence  $V_T$ ) or respiratory rate (RR). When  $\text{CO}_2$  rises, an increase in respiratory muscular effort usually precedes an increase in RR. Nevertheless, this depends also on a competent respiratory neuromuscular function as well as the lung mechanics.

While both are signs of respiratory distress, the RR and  $V_T$  themselves can furthermore induce/aggravate lung injury whether during MV (VILI) or spontaneous breathing (SILI). Furthermore, vigorous respiratory muscle contractions can precipitate myotrauma and ventilator induced diaphragmatic dysfunction (VIDD), increases the oxygen cost of breathing and consumes a significant part of the cardiac output. Altogether, these risk exhaustion and respiratory arrest. While the threshold to escalate support remains clinically based and difficult to define by simple parameters, undoubtedly a sound monitoring strategy can be of great help.




Figure 9: ( $P_{es}$ ) in  $\text{cmH}_2\text{O}$  on x-axis versus tidal volume in ml on y-axis. Green dashed line represents the chest wall compliance. Red shaded area is the Campbell diagram representing the inspiratory work of breathing. From reference 85.

### **15. Pendelluft phenomenon**

While ventilators display the  $V_t$  they deliver to the lungs, Pendelluft phenomenon refers to the exchange of volume in between the alveoli (movement of air within the lung from alveoli to alveoli without a change in VT). It usually happens in injured lungs during spontaneous effort from non-dependent to dependent areas. The reason is the non-uniform transmission of the diaphragmatic contraction (differential  $P_{pl}$ ) during spontaneous efforts. It follows that some alveoli will be subject to more strain, even with a VT within the target range. This highlights the value of regional monitoring of ventilation (see EIT in respiratory imaging).

### **16. Neuromuscular function**

A competent neuromuscular function is pivotal to counterbalance the respiratory load. It needs to be monitored probably from the start till the cure of ARF; nevertheless, 2 occasions are of greatest interest: (1) in the deteriorating patient leading to respiratory failure, and (2) when worsening or de novo muscle weakness risks weaning failure from MV. In the former, emergency situations herald effective monitoring during patients' distress and usually out of ICU. Following the institution of mechanical ventilation, diaphragmatic injury can happen after a few days of controlled ventilation which is associated with prolonged weaning and worse outcome.

Monitoring during weaning was much more investigated. The shift from controlled mode (total work done by the ventilator) to an assisted mode (e.g., pressure support ventilation, PSV) means the patient starts to share the load with the ventilator, which requires close observation. Both under and over assistance are associated with harm, and VILI/SILI are still possible as lungs usually did not heal completely. A mismatch between the patient and ventilator duration of insufflation can lead to double triggering or early cycling (PVA). So, the challenge here is to set the appropriate amount and duration of support, impossible without real time information about the respiratory drive, lung mechanics and most importantly the neuromuscular function.

#### **A. Airway Occlusion Pressure at 100 msec ( $P_{0.1}$ )**

$P_{0.1}$  refers to the  $P_{aw}$  at 100 msec after the start of inspiration while the respiratory circuit is occluded (Figure 9). The perception of any respiratory load is absent at that time point. As the airway is occluded, both the flow and volume equal zero, excluding any contribution of lung mechanics (i.e., resistance and elastance omitted). Its better correlation with the end tidal  $CO_2$  rather than the minute ventilation supports its correlation to the respiratory drive.  $P_{0.1}$  remains preserved even in muscle weakness if the patient can trigger. Last, while intrinsic PEEP can

theoretically cause a delay in triggering,  $P_{0.1}$  was validated in case of hyperinflation with good accuracy.

Most ventilators provide  $P_{0.1}$  monitoring ideally through end-expiratory occlusion despite its site (airway vs ventilatory circuit) and duration can differ affecting its value. Classically, it is performed as an intermittent manoeuvre and it is recommended to average 3-5 values. Some ventilators provide continuous breath per breath monitoring which permits the observation of the trend over time and with different levels of support. In these cases, occlusion is not possible with every breath and the value is extrapolated from the pressure curve. This risks inaccuracy (usually underestimation) and the operators should be mindful of the different methodologies between devices and their effect on the reference range.

$P_{0.1}$  ranges 0.5-1.5 cmH<sub>2</sub>O in healthy people, while values  $>3.5$  cmH<sub>2</sub>O was generally used to indicate increased drive during MV. It is of particular value during the sedation weaning and transition from controlled to assisted ventilation. A high  $P_{0.1}$  was shown to be associated with weaning failure, and in hyperinflated lungs, it can guide PEEP titration (drop in  $P_{0.1}$  indicates less intrinsic PEEP).  $P_{0.1}$  can help to differentiate double triggering (high value) from reverse triggering (low value).

### **B. Electrical activity of the diaphragm (EAdi)**

Monitoring the electrical activity of the respiratory neuromuscular function (electromyography) can be achieved via electrodes placed on the skin (surface electromyography, sEMG) or mounted on a nasogastric tube. While sEMG carries the additional benefit of monitoring the extra-diaphragmatic respiratory muscles, being non-invasive, but its application is technically challenging and needs high level of expertise limiting its use currently to research.

The diaphragm contributes 75-90% of the respiratory work. The electrical activity of the diaphragm (EAdi) is the sum of the diaphragmatic muscle fibers action potentials. The captured amplitude represents the diaphragmatic neural excitation (phrenic nerve activity) partially reflecting the central drive and its translation into a contraction magnitude per breath. The normal value in healthy people at rest is approximately 10  $\mu$ V.

The commercially available Neurally Adjusted Ventilatory Assist (NAVA) uses EAdi for better synchronization of the mechanical breath with the patient's neural inspiration and expiration.

EAdi is readily displayed on the ventilator screen and the operator sets a gain constant called the NAVA level which the ventilator translates into  $P_{aw}$ .

$$P_{aw} \text{ (cmH}_2\text{O)} = \text{EAdi (}\mu\text{V)} \cdot \text{NAVA level (cmH}_2\text{O}/\mu\text{V)}$$

Monitoring EAdi provides useful information to the clinician during weaning from MV, as it was shown that an early increase in its value was associated with weaning failure. The pressure generated by the respiratory muscles, usually measured by esophageal balloons, showed close correlation with EAdi in the single patient.

The relation between the electric activity and the generated pressure and volume can provide even more insight. Neuromechanical efficiency (NME) offers an estimation of the inspiratory effort and is calculated as  $\Delta P_{aw}$  divided by EAdi during an end-expiratory occlusion. Last, the ability to transform the electric activity into tidal volume is called Neuro-ventilatory efficiency index (NVE).

$$\text{NVE (ml}/\mu\text{V)} = \text{VT/EAdi}$$

EAdi/NAVA use can be affected by neurological injury of either the central respiratory centers or phrenic nerves. It measures only the diaphragmatic function, and liable to technical pitfalls including catheter positioning and signal interference (e.g., cardiac oscillations). It can be used during invasive or non-invasive ventilation and can help to synchronize the triggering and cycling of the breath.

### **C. Trans-diaphragmatic pressure (Pdi)**

Measuring the diaphragmatic force is another way to assess the diaphragmatic function. As Pressure equals force divided by the surface area, the pressure difference across the diaphragm is a valid measurement.

Transdiaphragmatic pressure (Pdi) is the difference between the pleural pressure and abdominal pressure

$$\mathbf{P_{di} = P_{pl} - P_{ab}}$$

An esophageal balloon with supplementary gastric balloon can measure the surrogate values of the  $P_{es}$  and gastric pressure ( $P_{ga}$ ):

$$\mathbf{P_{di} = P_{es} - P_{ga}}$$

Despite  $P_{di}$  remains a reference tool, its use is largely confined to research. It can be affected by the abdominal muscle contraction and lung volume (the change of the position and curvature of the diaphragm affect its length and hence its force).

#### **D. Respiratory muscle pressure ( $P_{mus}$ )**

Assessing the pressure generated by the respiratory muscles is the cornerstone to assess neuromuscular function.  $P_{aw}$  is readily shown on the ventilator screen but does not reflect  $P_{mus}$ . In fact, many derived measurements can be used bedside.

$P_{mus}$  is the reference for muscular assessment and can be computed from the motion equation if other variables are known (see above). This requires esophageal pressure monitoring to calculate the chest wall elastance ( $E_{cw}$ ) by dividing the  $\Delta P_{es}$  by the VT ( $E_{cw} = \Delta P_{es}/VT$ ). When the airway is occluded at the end of inspiration, the muscles become relaxed and their pressure is transmitted to the inner structures, increasing  $P_{pl}$ ,  $P_{alv}$  and  $P_{ao}$ . Then the  $P_{mus}$  can be calculated.

$$\mathbf{P_{mus} = V_t \times E_{cw} - \Delta P_{es}}$$

Such airway occlusion at the end of inspiration is also the basis to calculate the pressure muscle index (PMI) (see below).  $P_{mus}$  normally ranges between 3-15 cm H<sub>2</sub>O.

#### **E. Pressure Muscle Index (PMI):**

During PSV, the patient's muscles are generating a negative pressure usually measured by esophageal manometry. End-inspiratory hold was classically used to abolish the effect of flow and resistive pressure (static vs dynamic mechanics). However, it was shown that performing an end-inspiratory hold during PSV leading to muscle relaxation can reveal the pressure exerted by the muscles. To be valid, the maneuver needs to create a new constant plateau (flat plateau > 2 sec) with no visible inspiratory or expiratory effort. PMI can be measured as the difference between the end-inspiratory plateau and the sum of PS and PEEP. (Figure 10).



$$\mathbf{PMI = P_{end-insp\ plateau} - (PS + PEEP)}$$

Studies showed good correlation with  $P_{mus}$ , PTP, and other esophageal pressure-derived effort variables. The normal value ranges 0-2 cmH<sub>2</sub>O, while a value > 6 cmH<sub>2</sub>O denotes under-assistance and < 2 cmH<sub>2</sub>O means over-assistance.

#### **F. Maximum Inspiratory Pressure (MIP) and Maximum Expiratory Pressure (MEP)**

MIP ( $P_{I\max}$ ) and MEP ( $P_{E\max}$ ) are the peak pressures reached during maximum inspiration and expiration respectively. They need patient cooperation for valid testing and are recommended by respiratory societies as a global marker of muscular strength.

MIP indicates inspiratory muscle strength. In ventilated patients, end expiratory occlusion is used to measure the generated inspiratory pressure and hence the name: Occlusion pressure ( $P_{occ}$ ) (see below).

MEP also indicates muscle strength but it is particularly important to assess coughing and secretion clearance, a fundamental prerequisite for successful extubation especially in patients with muscle weakness. The classical MEP (cooperative sitting patient with a nose clip) is not possible in the ventilated ICU patient. Instead, it was suggested to measure the pressure and flow during induced cough (Cough peak flow (CPF)). CPF value depends on both the inspiratory and expiratory muscle function, lung and airway mechanics, age and the glottic closure (impossible in intubated patients). This is usually tested after a successful spontaneous breathing trial (SBT) in non-sedated patients when the best flow velocity can be recorded by the ventilator (or by an external device). A value >60 L/min at the ETT level was suggested to predict successful extubation. It can also be used before decannulation of tracheostomies.

#### **G. Occlusion Pressure ( $P_{occ}$ )**

In the ICU, many ventilators offer the option to perform the maximum inspiratory pressure (MIP) in intubated patients, but it requires cooperation and repeated measurements. It is another occlusion pressure (like  $P_{0.1}$ ) requiring end-expiratory hold. While  $P_{0.1}$  captures the pressure earlier at 0.1 sec,  $P_{occ}$  will represent the delta change in pressure between end expiration (PEEP) and the negative peak (trough) during occlusion. (Figure 10)

$$\mathbf{P_{occ} = P_{trough} - PEEP}$$

Technical differences between manufacturers do exist which includes the duration of the occlusion and the starting volume (residual volume vs functional residual capacity). The normal range is considered between -3 and -15 cmH<sub>2</sub>O. A high value suggest weakness is less likely, while a low value can be biased by technical factors or poor effort.

It represents the muscular effort and total stress of the lungs but the main interest in P<sub>occ</sub> may stem from the validated calculation of P<sub>mus</sub> from its value:

$$P_{mus} = K \cdot P_{occ} = -3/4 \cdot P_{occ}$$

K is a constant which was shown to be equal to  $-3/4$

Furthermore, knowing P<sub>mus</sub> means that the dynamic TPP can be calculated:

$$TPP_{(dyn)} = (P_{peak} - PEEP) - (2/3 P_{occ})$$

Online calculator can be found at: <https://societymechanicalventilation.org/calculators/>




Figure10: Ventilator interface showing inspiratory effort. Key features are labelled for measurement of PMI (during inspiratory hold), P<sub>occ</sub> and P<sub>0.1</sub> (measuring during expiratory hold). Arrows denote the beginning and end of the maneuver.

## **Conclusion**

A personalized approach of respiratory support is not possible without a multifaceted monitoring strategy. Different tools can be used concomitantly or subsequently in different phases of the disease, particularly in the most complex cases, namely (but not exclusively) ARDS.

The generated data should be interpreted not only correctly but also in a timely manner. To this end, a deep understanding of the respiratory physiology, underlying pathology, monitoring principles and the limits of a growing number of tools is becoming a challenging mission for clinicians. The inclination to use machine learning and artificial intelligence for automatic closed loop adjustment is promising but will not replace human supervision at least in the near future.

## References

1. GBD 2019 Chronic Respiratory Diseases Collaborators. Global burden of chronic respiratory diseases and risk factors, 1990-2019: an update from the Global Burden of Disease Study 2019. *EClinicalMedicine* 2023; 59:101936.
2. GBD 2021 Upper Respiratory Infections Otitis Media Collaborators. Global, regional, and national burden of upper respiratory infections and otitis media, 1990-2021: a systematic analysis from the Global Burden of Disease Study 2021. *Lancet Infect Dis* 2025; 25(1):36-51.
3. Vitazkova D, Foltan E, Kosnacova H, et al. Advances in respiratory monitoring: A comprehensive review of wearable and remote technologies. *Biosensors (Basel)* 2024; 14(2):90.
4. Roshdy A. Respiratory monitoring during mechanical ventilation: The present and the future. *J Intensive Care Med* 2023; 38(5):407-417.
5. Brochard L, Martin GS, Blanch L, et al. Clinical review: Respiratory monitoring in the ICU - a consensus of 16. *Crit Care* 2012; 16(2):219.
6. Daoud EG, Katigbak R, Ottochian M. Accuracy of the ventilator automated displayed respiratory mechanics in passive and active breathing conditions: A bench study. *Respir Care* 2019; 64(12):1555-1560.
7. Bellani G, Laffey JG, Pham T, et al; LUNG SAFE Investigators; ESICM Trials Group. Epidemiology, patterns of care, and mortality for patients with acute respiratory distress syndrome in intensive care units in 50 countries. *JAMA*. 2016; 315(8):788-800.
8. Coppola S, Pozzi T, Catozzi G, et al. Clinical performance of Spo<sub>2</sub>/Fio<sub>2</sub> and Pao<sub>2</sub>/Fio<sub>2</sub> ratio in mechanically ventilated acute respiratory distress syndrome patients: A retrospective study. *Crit Care Med*. Published online March 3, 2025.
9. Junhai Z, Jing Y, Beibei C, et al. The value of ROX index in predicting the outcome of high flow nasal cannula: a systematic review and meta-analysis. *Respir Res* 2022; 23(1):33.
10. Roca O, Caralt B, Messika J, et al. An index combining respiratory rate and oxygenation to predict outcome of nasal high-flow therapy. *Am J Respir Crit Care Med* 2019; 199(11):1368–1376.
11. Roshdy A. Echodynamics: interpretation, limitations, and clinical integration!. *J Intensive Care Med* 2018; 33(8):439-446.
12. Kollisch-Singule M, Satalin J, Blair SJ, et al. Mechanical ventilation lessons learned from alveolar micromechanics. *Front Physiol* 2020; 11:233.
13. Albert SP, DiRocco J, Allen GB, et al. The role of time and pressure on alveolar recruitment. *J Appl Physiol (1985)* 2009; 106(3):757-765.

14. Creamer KM, McCloud LL, Fisher LE, et al. Closing pressure rather than opening pressure determines optimal positive end-expiratory pressure and avoids overdistention. *CHEST* 1999; 116(1 Suppl):26S-27S.
15. Katayama S, Iotti GA, Tonai K. Airway opening pressure in an obese ARDS patient: just the start of airway recruitment. *J Anesth Analg Crit Care* 2026; 3(1):13.
16. Chen L, Del Sorbo L, Grieco DL, et al. Airway closure in acute respiratory distress syndrome: an underestimated and misinterpreted phenomenon. *Am J Respir Crit Care Med* 2018; 197(1):132-136.
17. Crotti S, Mascheroni D, Caironi P, et al. Recruitment and derecruitment during acute respiratory failure a clinical study. *Am J Respir Crit Care Med*. 2001;164:131–140.
18. <https://www.getinge.com/int/company/news/press-releases/2025/5009600-Getinge-introduces-Neural-Pressure-Support-for-the-Servo-u/> as accessed on the 19.03.2025
19. Carrera M, Urrutia JG, Ardariz CB, et al. Maximal expiratory pressure compared with maximal expiratory pressure during induced cough as a predictor of extubation failure. *Crit Care Sci* 2023; 35(1):37-43.
20. Brennan M, McDonnell MJ, Duignan N, et al. The use of cough peak flow in the assessment of respiratory function in clinical practice- A narrative literature review. *Respir Med* 2022; 193:106740.
21. Noeman-Ahmed Y, Gokaraju S, Powrie DJ, et al. Predictors of CPAP outcome in hospitalized COVID-19 patients. *Respirology* 2020; 25(12):1316-1319.
22. Matthay MA, Arabi Y, Arroliga AC, et al. A New global definition of acute respiratory distress syndrome. *Am J Respir Crit Care Med* 2024; 209(1):37-47.
23. De Santis V, Singer M. Tissue oxygen tension monitoring of organ perfusion: rationale, methodologies, and literature review. *Br J Anaesth* 2015; 115(3):357-365.
24. Figaji AA, Zwane E, Graham Fieggen A, et al. The effect of increased inspired fraction of oxygen on brain tissue oxygen tension in children with severe traumatic brain injury. *Neurocrit Care* 2010; 12(3):430-437.
25. Teboul JL, Monnet X, De Backer D. Should We Abandon Measuring SvO<sub>2</sub> or ScvO<sub>2</sub> in Patients with Sepsis? In: Vincent, JL. (eds) *Annual Update in Intensive Care and Emergency Medicine* 2019. *Annual Update in Intensive Care and Emergency Medicine*. Springer, Cham.
26. Chetana Shanmukhappa S, Lokeshwaran S. Venous Oxygen Saturation. In: *StatPearls*. Treasure Island (FL): StatPearls Publishing; September 10, 2024.

27. Sakuraya M, Douno E, Iwata W, et al. Accuracy evaluation of mainstream and sidestream end-tidal carbon dioxide monitoring during noninvasive ventilation: a randomized crossover trial (MASCAT-NIV trial). *J Intensive Care* 2022; 10(1):17.
28. Lecompte-Osorio P, Pearson SD, Pieroni CH, et al. Bedside estimates of dead space using end-tidal CO<sub>2</sub> are independently associated with mortality in ARDS. *Crit Care* 2021;25(1):333.
29. Daoud EG, Franck CL. Alveolar mechanics: A new concept in respiratory monitoring. *J Mech Vent* 2022; 3(4):178-188.
30. Akoumianaki E, Maggiore SM, Valenza F, et al. The application of esophageal pressure measurement in patients with respiratory failure. *Am J Respir Crit Care Med* 2014; 189:520–531.
31. Cabello B, Mancebo J. Work of breathing. *Intensive Care Med* 2006; 32:1311–1314.
32. Apigo M, Schechtman J, Dhliwayo N, et al. Development of a work of breathing scale and monitoring need of intubation in COVID-19 pneumonia. *Crit Care* 2020; 24(1):477.
33. J Brady Scott, Ramandeep Kaur. Monitoring breathing frequency, pattern, and effort. *Respire Care* 2020; 65(6):793-806.
34. Gattinoni L, Giosa L, Bonifazi M, et al. Targeting transpulmonary pressure to prevent ventilator-induced lung injury. *Expert Rev Respir Med* 2019; 13(8):737–746.
35. <https://www.hamilton-medical.com/de/Article-page~knowledge-base~3092b082-f32b-4106-9720-ce4d94a308b1~Can-I-measure-work-of-breathing-for-patient-and-ventilator~.html> as accessed 24 January 2025.
36. Ince C. Personalized physiological medicine. *Crit Care* 2017; 21(Suppl 3):308.
37. Fritsch SJ, Cecconi M. Setting the ventilator with AI support: challenges and perspectives. *Intensive Care Med*. Published online January 20, 2025.
38. Harms CA, Wetter TJ, McClaran SR, et al. Effects of respiratory muscle work on cardiac output and its distribution during maximal exercise. *J Appl Physiol* (1985) 1998 ;85(2):609-618.
39. Goligher EC, Jonkman AH, Dianti J, et al. Clinical strategies for implementing lung and diaphragm-protective ventilation: avoiding insufficient and excessive effort. *Intensive Care Med* 2020; 46(12):2314-2326.
40. Shi ZH, Jonkman A, de Vries H, et al. Expiratory muscle dysfunction in critically ill patients: towards improved understanding. *Intensive Care Med* 2019; 45(8):1061-1071.
41. Goligher EC, Fan E, Herridge MS, et al. Evolution of diaphragm thickness during mechanical ventilation, impact of inspiratory effort. *Am J Respir Crit Care Med* 2015; 192:1080–1088.

42. Grasselli G, Calfee CS, Camporota L, et al. ESICM guidelines on acute respiratory distress syndrome: definition, phenotyping and respiratory support strategies. *Intensive Care Med* 2023; 49:727–759.
43. Carreaux G, Millán-Guilarte T, De Prost N, et al. Failure of noninvasive ventilation for de novo acute hypoxemic respiratory failure: Role of tidal volume. *Crit Care Med* 2016 ;44(2):282-290.
44. Akoumianaki E, Vaporidi K, Georgopoulos D. The injurious effects of elevated or nonelevated respiratory rate during mechanical ventilation. *Am J Respir Crit Care Med* 2019; 199(2):149-157.
45. Steven Cruickshank, Nicola Hirschauer, The alveolar gas equation. *Continuing Education in Anaesthesia Critical Care & Pain* 2004; 1: 24–27
46. Buist M, Bernard S, Nguyen TV, et al. Association between clinically abnormal observations and subsequent in-hospital mortality: a prospective study. *Resuscitation* 2004; 62:137–141.
47. Cretikos MA, Bellomo R, Hillman K, et al. Respiratory rate: the neglected vital sign. *Med J Aust* 2008; 188:657–659.
48. Fabry B. How to minimize mechanical power during controlled mechanical ventilation. *Intensive Care Med Exp* 2024; 12(1):114.
49. Webb HH, Tierney DF. Experimental pulmonary edema due to intermittent positive pressure ventilation with high inflation pressures. Protection by positive end-expiratory pressure. *Am Rev Respir Dis* 1974; 110:556–565.
50. Dreyfuss D, Soler P, Basset G, et al. High inflation pressure pulmonary edema. Respective effects of high airway pressure, high tidal volume, and positive end-expiratory pressure. *Am Rev Respir Dis* 1988; 137:1159–1164.
51. Singer M, Young PJ, Laffey JG. et al. Dangers of hyperoxia. *Crit Care* 2021; 25(1):440.
52. Buell KG, Spicer AB, Casey JD, et al. Individualized treatment effects of oxygen targets in mechanically ventilated critically ill adults. *JAMA* 2024; 331(14):1195–1204.
53. Semler MW, Casey JD, Lloyd BD, et al. Oxygen-saturation targets for critically ill adults receiving mechanical ventilation. *N Engl J Med* 2022; 387(19):1759-1769.
54. Martin DS, Grocott MPW. Heterogeneity of treatment effect: the case for individualising oxygen therapy in critically ill patients. *Crit Care* 2025; 29(1):50.
55. Cressoni M, Chiumello D, Carlesso E, et al. Compressive forces and computed tomography–derived positive end expiratory pressure in acute respiratory distress syndrome. *Anesthesiology* 2014; 121:572–581.
56. Yoshida T, Amato MBP, Grieco DL, et al. Esophageal manometry and regional transpulmonary pressure in lung injury. *Am J Respir Crit Care Med* 2018; 197(8):1018-1026.

57. Loring SH, Topulos GP, Hubmayr RD. Transpulmonary pressure: the importance of precise definitions and limiting assumptions. *Am J Respir Crit Care Med* 2016; 194(12):1452-1457.
58. Shokry M, Yamasaki K, Daoud EG. Can you calculate the total respiratory, lung, and chest wall respiratory mechanics? *J Mech Vent* 2020; 1(1):24-25.
59. Tonetti T, Vasques F, Rapetti F, et al. Driving pressure and mechanical power: new targets for VILI prevention. *Ann Transl Med* 2017; 5(14): 286.
60. Amato MB, Meade MO, Slutsky AS, et al. Driving pressure and survival in the acute respiratory distress syndrome. *N Engl J Med* 2015; 372:747-755.
61. Depta F, Chiofolo CM, Chbat NW, et al. Six methods to determine expiratory time constants in mechanically ventilated patients: a prospective observational physiology study. *Intensive Care Med Exp* 2024; 12(1):25.
62. Giosa L, Collins PD, Shetty S, et al. Bedside assessment of the respiratory system during invasive mechanical ventilation. *J Clin Med* 2024; 13(23):7456.
63. Henderson WR, Dominelli PB, Molgat-Seon Y, et al. Effect of tidal volume and positive end-expiratory pressure on expiratory time constants in experimental lung injury. *Physiol Rep* 2016; 4(5):e12737.
64. Hamilton-medical. Using ventilator waveform analysis to detect patient-ventilator asynchronies. The reference card. [https://www.hamilton-medical.com/en\\_US/Landing-pages/Patient-ventilator-asynchrony-cheatsheet.html](https://www.hamilton-medical.com/en_US/Landing-pages/Patient-ventilator-asynchrony-cheatsheet.html) (accessed 12 Feb 2025).
65. de Haro C, Ochagavia A, López-Aguilar J, et al. Patient-ventilator asynchronies during mechanical ventilation: current knowledge and research priorities. *Intensive Care Med Exp* 2019; 7(Suppl 1):43.
66. Thille AW, Rodriguez P, Cabello B, et al. Patient-ventilator asynchrony during assisted mechanical ventilation. *Intensive Care Med* 2006; 32(10):1515-1522.
67. Beitler JR, Sands SA, Loring SH, et al. Quantifying unintended exposure to high tidal volumes from breath stacking dyssynchrony in ARDS: the BREATHE criteria. *Intensive Care Med* 2016; 42(9):1427-1436.
68. Gattinoni L, Mascheroni D, Torresin A, et al. Morphological response to positive end expiratory pressure in acute respiratory failure. Computerized tomography study. *Intensive Care Med* 1986; 12(3):137-142.
69. Gattinoni L, Marini JJ, Pesenti A. et al. The "baby lung" became an adult. *Intensive Care Med* 2016; 42:663–673. (2016).



70. Writing Group for the Alveolar Recruitment for Acute Respiratory Distress Syndrome Trial (ART) Investigators, Cavalcanti AB, Suzumura ÉA, et al. Effect of lung recruitment and titrated positive end-expiratory pressure (PEEP) vs low PEEP on mortality in patients with acute respiratory distress syndrome: A randomized clinical trial. *JAMA* 2017; 318(14):1335-1345.
71. Giovanazzi S, Nocera D, Catozzi G, et al. Assessment of recruitment from CT to the bedside: challenges and future directions. *Crit Care* 2025; 29(1):64.
72. Chen L, Del Sorbo L, Grieco DL, et al. Potential for lung recruitment estimated by the recruitment-to inflation ratio in acute respiratory distress syndrome a clinical trial. *Am J Respir Crit Care Med* 2020; 201:178–187.
73. Chiumello D, Marino A, Brioni M, Cigada I, Menga F, Colombo A, et al. Lung recruitment assessed by respiratory mechanics and computed tomography in patients with acute respiratory distress syndrome what is the relationship? *Am J Respir Crit Care Med* 2016; 193:1254–1263.
74. Murgolo F, Grieco DL, Spadaro S, et al. Recruitment to-inflation ratio reflects the impact of peep on dynamic lung strain in a highly recruitable model of ARDS. *Ann Intensive Care* 2024; 14:106.
75. Cour M, Biscarrat C, Stevic N, et al. Recruitment-to-inflation ratio measured with modern intensive care unit ventilators: How accurate is it? *Crit Care* 2022; 26(1):85.
76. Del Sorbo L, Manuel Tisminetzky Lu, Chen LB, et al. Association of lung recruitment and change in recruitment-to inflation ratio from supine to prone position in acute respiratory distress syndrome. *Crit Care* 2023; 27(1):140.
77. Escolar JD, Escolar A. Lung hysteresis: a morphological view. *Histol Histopathol* 2004; 19(1):159-166.
78. Protti A, Andreis DT, Milesi M, et al. Role of amplitude and rate of deformation in ventilator-induced lung injury. *Intensive Care Med Exp* 2014; 2(Suppl 1):P75.
79. Marini JJ. Evolving concepts for safer ventilation. *Crit Care* 2019; 23(Suppl 1):114.
80. [https://www.hamilton-medical.com/en\\_GB/Products/Technologies/P-V-Tool.html](https://www.hamilton-medical.com/en_GB/Products/Technologies/P-V-Tool.html). Accessed 14 Feb 2025
81. Demory D, Arnal JM, Wysocki M, et al. Recruitability of the lung estimated by the pressure volume curve hysteresis in ARDS patients. *Intensive Care Med* 2008; 34(11):2019-2025.
82. Chiumello D, Arnal JM, Umbrello M, et al. Hysteresis and lung recruitment in acute respiratory distress syndrome patients: a CT scan study. *Crit Care Med* 2020; 48:1494–1502.
83. Nakayama R, Bunya N, Katayama S, et al. Correlation between the hysteresis of the pressure–volume curve and the recruitment-to-inflation ratio in patients with coronavirus disease 2019. *Ann. Intensive Care* 2022; 12(1):106.

84. Sundaram M, Karthika M. Respiratory mechanics: To balance the mechanical breaths!!. *Indian J Crit Care Med* 2021; 25(1):10-11.
85. Shokry M, Simonpietri M, Yamasaki K, Daoud EG. Calculating the work of breathing during mechanical ventilation. *J Mech Vent* 2021; 2(2):71-72.
86. Jimenez JV, Munroe E, Weirauch AJ, et al. Electric impedance tomography-guided PEEP titration reduces mechanical power in ARDS: a randomized crossover pilot trial. *Crit Care* 2023; 27(1):21.
87. Bouhemad B, Brisson H, Le-Guen M, et al. Bedside ultrasound assessment of positive end-expiratory pressure-induced lung recruitment. *Am J Respir Crit Care Med* 2011; 183(3):341-347.
88. Tusman G, Acosta CM, Costantini M. Ultrasonography for the assessment of lung recruitment maneuvers. *Crit Ultrasound J* 2016; (1):8.
89. Tusman G, Acosta CM, Nicola M, et al. Real-time images of tidal recruitment using lung ultrasound. *Crit Ultrasound* 2015; 7(1):19.
90. Mansouri S, Alharbi Y, Haddad F, et al. Electrical Impedance Tomography - Recent applications and developments. *J Electr Bioimpedance* 2021; 12(1):50-62.
91. Scaramuzzo G, Pavlovsky B, Adler A, et al. Electrical impedance tomography monitoring in adult ICU patients: state-of-the-art, recommendations for standardized acquisition, processing, and clinical use, and future directions. *Crit Care* 2024; 28:377.
92. Xu M, He H, Long Y. Lung perfusion assessment by bedside electrical impedance tomography in critically ill patients. *Front Physiol* 2021; 12:748724.
93. Elwan MH, Roshdy A, Elsharkawy EM, et al. Can passive leg raise predict the response to fluid resuscitation in ED? *BMC Emerg Med* 2022; 22(1):172.
94. Yoshida T, Piraino T, Lima CAS, et al. Regional ventilation displayed by electrical impedance tomography as an incentive to decrease positive end-expiratory pressure. *Am J Respir Crit Care Med* 2019; 200(7):933-937.
95. Costa EL, Borges JB, Melo A, et al. Bedside estimation of recruitable alveolar collapse and hyperdistension by electrical impedance tomography. *Intensive Care Med* 2009; 35(6):1132-1137.
96. Wang Yx, Zhong M, Dong Mh, et al. Prone positioning improves ventilation–perfusion matching assessed by electrical impedance tomography in patients with ARDS: a prospective physiological study. *Crit Care* 2022; 26(1):154.
97. Franchineau G, Bréchet N, Hekimian G, et al. Prone positioning monitored by electrical impedance tomography in patients with severe acute respiratory distress syndrome on venovenous ECMO. *Ann Intensive Care* 2020; 10(1):12.

98. Zhou T, Zhang HP, Chen WW, et al. Cuff-leak test for predicting postextubation airway complications: a systematic review. *Journal of Evidence-Based Medicine*, 2011; 4:242-254.
99. Sahetya SK, Goligher EC, Brower RG. Fifty years of research in ards. setting positive end-expiratory pressure in acute respiratory distress syndrome. *Am J Respir Crit Care Med* 2017; 195(11):1429-1438.
100. Web HH, Tierney DF. Experimental pulmonary edema due to intermittent positive pressure ventilation with high inflation pressures. protection by positive end-expiratory pressure. *Am Rev Respir Dis* 1974; 110(5):556-565.
101. Zersen KM. Setting the optimal positive end-expiratory pressure: a narrative review. *Front Vet Sci* 2023; 10:1083290.
102. Bellani G, Laffey JG, Pham T, et al. Epidemiology, patterns of care, and mortality for patients with acute respiratory distress syndrome in intensive care units in 50 countries. *JAMA* 2016; 315:788–800.
100. Ball, L., Talmor, D. & Pelosi, P. Transpulmonary pressure monitoring in critically ill patients: pros and cons. *Crit Care* 2024; 28(1):177.
101. Edginton S, Kruger N, Stelfox HT, et al. Methods for determining optimal positive end-expiratory pressure in patients undergoing invasive mechanical ventilation: a scoping review. *Can J Anaesth*. 2024; 71(11):1535-1555.
102. Shono A, Kotani T. Clinical implication of monitoring regional ventilation using electrical impedance tomography. *J Intensive care* 2019; 7:4.
103. Muders T, Luepschen H, Zinserling J, et al. Tidal recruitment assessed by electrical impedance tomography and computed tomography in a porcine model of lung injury. *Crit Care Med* 2012; 40(3):903–911.
104. Tomacic V, Cornejo R. Lung monitoring with electrical impedance tomography: technical considerations and clinical applications. *J Thorac Dis* 2019; 11(7):3122-3135.
105. Costa EL, Gonzalez R, Amato MB. Electrical impedance tomography. *Curr Opin Crit Care* 2009; 15:18-24.
106. Anne-Kathrin Brill. How to avoid interface problems in acute noninvasive ventilation. *Breathe* 2014; 10(3):230-242.
107. Chen Q, Yu X, Chen Y, et al. Intermittent versus continuous control of tracheal cuff pressure in patients undergoing mechanical ventilation: A systematic review and meta-analysis. *J Clin Nurs* 2023; 32(15-16):4283-4294.

108. Amato MB, Meade MO, Slutsky AS, et al. Driving pressure and survival in the acute respiratory distress syndrome. *N Engl J Med* 2015; 372(8):747-755.
109. Zaidi SF, Shaikh A, Khan DA, et al. Driving pressure in mechanical ventilation: A review. *World J Crit Care Med* 2024; 13(1):88385.
110. Elliott MW. The interface: crucial for successful noninvasive ventilation. *Eur Respir J* 2004; 23(1):7-8.
111. Schwartz, Evan A, Bronshteyn YS, et al. Ventilator Stress Index: An Intensive Care Unit Tool That Anesthesiologists Should Know. *J Cardiothorac Vasc Anesth.* 2024 Dec;38(12):3252-3256.
112. Gattinoni L, Collino F, Camporota L. Mechanical power: meaning, uses and limitations. *Intensive Care Med* 2023; 49(4):465-467.
113. Huhle R, Serpa-Neto A, Schultz MJ, et al. Is mechanical power the final word on ventilator-induced lung injury? No. *Ann Transl Med* 2018; 6(19):394.
114. Acute Respiratory Distress Syndrome Network. Ventilation with lower tidal volumes as compared with traditional tidal volumes for acute lung injury and the acute respiratory distress syndrome. *N Engl J Med* 2000 ;342(18):1301-1308.
115. von Düring S, Parhar KKS, Adhikari NKJ, et al. Understanding ventilator-induced lung injury: The role of mechanical power. *J Crit Care* 2025; 85:154902.
116. Guerin C, Papazian L, Reignier J, et al. Effect of driving pressure on mortality in ARDS patients during lung protective mechanical ventilation in two randomized controlled trials. *Crit Care* 2016; 20:384.
117. Serpa Neto A, Deliberato RO, Johnson AEW, et al. Mechanical power of ventilation is associated with mortality in critically ill patients: an analysis of patients in two observational cohorts. *Intensive Care Med* 2018; 44(11):1914-1922.
118. Justin Seemueller. How to estimate mechanical power in volume and pressure control ventilation. [https://www.hamilton-medical.com/en\\_GB/News-Events/News-Archive/Article-page~knowledge-base~f40239e0-e478-43ad-ba33-121f81fe632b~How-to-estimate-mechanical-power-in-volume--and-pressure-control-ventilation~.html](https://www.hamilton-medical.com/en_GB/News-Events/News-Archive/Article-page~knowledge-base~f40239e0-e478-43ad-ba33-121f81fe632b~How-to-estimate-mechanical-power-in-volume--and-pressure-control-ventilation~.html). Accessed 20 Feb 2025
119. Takaoka K, Toma S, Lee P, et al. A comparative analysis of mechanical power and Its components in pressure-controlled ventilation mode and AVM-2 mode. *J Mech Vent* 2023; 4(4):130-140.
120. Yoshida T, Torsani V, Gomes S, et al. Spontaneous effort causes occult pendelluft during mechanical ventilation. *Am J Respir Crit Care Med* 2013; 188(12):1420-1427.
121. Jonkman, AH, de Vries HJ, Heunks LMA. Physiology of the respiratory drive in icu patients: implications for diagnosis and treatment. *Crit Care* 2020; ;28(1):94.

122. Yoshida T, Amato MBP, Kavanagh BP. Understanding spontaneous vs. ventilator breaths: impact and monitoring. *Intensive Care Med* 2018; 44(12):2235–2238.
123. Georgopoulos D, Bolaki M, Stamatopoulou V, et al. Respiratory drive: a journey from health to disease. *J Intensive Care* 2024; 12:15.
124. Peñuelas O, Keough E, López-Rodríguez L, et al. Ventilator-induced diaphragm dysfunction: translational mechanisms lead to therapeutical alternatives in the critically ill. *Intensive Care Med Exp* 2019; 7(Suppl 1):48.
125. Vaporidi K, Akoumianaki E, Telias I, et al. Respiratory drive in critically ill patients: pathophysiology and clinical implications. *Am J Respir Crit Care Med* 2020; 201(1):20–32.
126. Whitelaw WA, Derenne JP, Milic-Emili J. Occlusion pressure as a measure of respiratory center output in conscious man. *Respir Physiol* 1975; 23:181–199.
127. Telias I, Damiani F, Brochard L. The airway occlusion pressure (P0.1) to monitor respiratory drive during mechanical ventilation: increasing awareness of a not-so-new problem. *Intensive Care Med* 2018; 44(9):1532-1535.
128. Mancebo J, Albaladejo P, Touchard D, et al. Airway occlusion pressure to titrate positive end-expiratory pressure in patients with dynamic hyperinflation. *Anesthesiology* 2000; 93:81–90.
129. Katayama S, Tonai K, Nunomiya S. Bias and precision of continuous P0.1 measurement by various ventilators: A simulation study. *Respir Care* 2023; 68(10):1393-1399.
130. Beloncle F, Piquilloud L, Olivier PY, et al. Accuracy of P0.1 measurements performed by ICU ventilators: a bench study. *Ann Intensive Care* 2019; 9:104.
131. Sklar MC, Burns K, Rittayamai N, et al. Effort to breathe with various spontaneous breathing trial techniques. A physiologic meta-analysis. *Am J Respir Crit Care Med* 2017; 195(11):1477–1485.
132. Rittayamai N, Beloncle F, Goligher EC, et al. Effect of inspiratory synchronization during pressure-controlled ventilation on lung distension and inspiratory effort. *Ann Intensive Care* 2017; 7:100.
133. Conti G, Cinnella G, Barboni E, et al. Estimation of occlusion pressure during assisted ventilation in patients with intrinsic PEEP. *Am J Respir Crit Care Med* 1996; 154(4 Pt 1):907-912.
134. Sant'ambrogio G, Camporesi E. Contribution of various inspiratory muscles to ventilation and the immediate and distant effect of diaphragmatic paralysis. *Acta Neurobiol Exp (Wars)* 1973; 33(1):401-409.
135. Diniz-Silva F, Miethke-Morais A, Alencar AM, et al. Monitoring the electric activity of the diaphragm during noninvasive positive pressure ventilation: a case report. *BMC Pulm Med* 2017; 17:91.

136. Jonkman AH, Warnaar RSP, Baccinelli W, et al. Analysis and applications of respiratory surface EMG: report of a round table meeting. *Crit Care* 2024; 28(1):2.
137. Navalesi P, Colombo D, Della Corte F. NAVA ventilation. *Minerva Anestesiol* 2010; 76(5):346-352.
138. Barwing J, Pedroni C, Olgemöller U, et al. Electrical activity of the diaphragm (EAdi) as a monitoring parameter in difficult weaning from respirator: a pilot study. *Crit Care* 2013; 17(4):R182.
139. Yao L, Zhou Z, Yu T, et al. The predictive value of neurally adjusted ventilatory assist indexes for the prognosis of patients with severe cerebral hemorrhage. *Eur J Med Res* 2024; 29:7.
140. Bellani G, Mauri T, Coppadoro A, et al. Estimation of patient's inspiratory effort from the electrical activity of the diaphragm. *Crit Care Med* 2013; 41(6):1483-1491.
141. Inata Y, Takeuchi M. Ventilator auto-triggering by cardiac electrical activity during noninvasive ventilation with neurally adjusted ventilatory assist. *Clin Case Rep* 2018; 6(7):1379-1380.
142. Brochard L. Transdiaphragmatic Pressure. In: Benito, S., Net, A. (eds) *Pulmonary Function in Mechanically Ventilated Patients. Update in Intensive Care and Emergency Medicine*, 1991; vol 13. Springer, Berlin, Heidelberg.
143. Umbrello M, Formenti P, Lusardi AC, et al. Oesophageal pressure and respiratory muscle ultrasonographic measurements indicate inspiratory effort during pressure support ventilation. *Br J Anaesth* 2020; 125(1):e148-e157.
144. Foti G, Cereda M, Banfi G, et al. End-inspiratory airway occlusion: a method to assess the pressure developed by inspiratory muscles in patients with acute lung injury undergoing pressure support. *Am J Respir Crit Care Med* 1997; 156:1210-1216.
145. Pavez N, Damiani LF. Inspiratory and expiratory pause during pressure support ventilation: Maneuvers that we should incorporate into clinical practice. *Med Intensiva (Engl Ed)*. 2022; 46(4):213-216.
146. Mauri T, Guérin C, Hubmayr R. The ten pressures of the respiratory system during assisted breathing. *Intensive Care Med* 2017; 43(10):1504-1506.
147. Mauri T, Cambiaghi B, Spinelli E, et al. Spontaneous breathing: a double-edged sword to handle with care. *Ann Transl Med* 2017; 5(14):292.
148. Natalini G, Buizza B, Granato A, et al. Non-invasive assessment of respiratory muscle activity during pressure support ventilation: accuracy of end-inspiration occlusion and least square fitting methods. *J Clin Monit Comput* 2021; 35:913–921.

149. Gao R, Yang Y, Zhang L, et al. Use of pressure muscle index to guide pressure support ventilation setting: a study protocol and statistical plan for a prospective randomised controlled proof-of-concept trial. *BMJ Open* 2024; 14:e082395.
150. Pmus Index - Plateau Pressure during pressure support ventilation. <https://rtmaven.com/pmi> as accessed on 23 Feb 2023
151. Bellani G, Grassi A, Sosio S et al. Plateau and driving pressure in the presence of spontaneous breathing. *Intensive Care Med* 2019; 45:97–98.
152. Goligher EC, Dres M, Fan E, et al. Mechanical Ventilation-induced Diaphragm Atrophy Strongly Impacts Clinical Outcomes. *Am J Respir Crit Care Med* 2018; 197(2):204-213.
153. American Thoracic Society/European Respiratory Society. ATS/ERS Statement on respiratory muscle testing. *Am J Respir Crit Care Med* 2002; 166(4):518-624.
154. Laveneziana P, Albuquerque A, Aliverti A, et al. ERS statement on respiratory muscle testing at rest and during exercise. *Eur Respir J* 2019; 53(6):1801214.
155. Bertoni M, Telias I, Urner M, et al. A novel non-invasive method to detect excessively high respiratory effort and dynamic transpulmonary driving pressure during mechanical ventilation. *Crit Care* 2019; 23:346.
156. Jaber S, Petrof BJ, Jung B, et al. Rapidly progressive diaphragmatic weakness and injury during mechanical ventilation in humans. *Am J Respir Crit Care Med* 2011; 183:364–371.
157. Schellekens WJM, van Hees HWH, Doorduyn J, et al. Strategies to optimize respiratory muscle function in ICU patients. *Crit Care* 2016; 20:103. (2016).
158. van Oosten JP, Akoumianaki E, Jonkman AH. Monitoring respiratory muscles effort during mechanical ventilation. *Curr Opin Crit Care* 2025; 31(1):12-20.
159. Hermans G, Demoule A, Heunks L. How I perform diaphragmatic ultrasound in the intensive care unit. *Intensive Care Med* 2024; 50:2175–2178.
160. Tuinman PR, Jonkman AH, Dres M, et al. Respiratory muscle ultrasonography: methodology, basic and advanced principles and clinical applications in ICU and ED patients—a narrative review. *Intensive Care Med* 2020; 46(4):594–605.
161. Cardenas LZ, Santana PV, Caruso P, et al. Diaphragmatic ultrasound correlates with inspiratory muscle strength and pulmonary function in healthy subjects. *Ultrasound Med Biol* 2018; 44:786–793.
162. Hu J, Hasan O, Shiraishi K, Hirao Y, Daoud EG. Comparison of estimation of inspiratory muscle effort using three common indices in various respiratory models, a bench study. *J Mech Vent* 2024; 5(4):119-125.